# OEC UroView<sup>®</sup> 2800

## Table Motion Troubleshooting Procedure

00-885283-01

*Advanced Service Information  
Company Confidential*

  
Contents

  
Installation

  
Service

  
Schematics

  
Periodic Maintenance

  
Illustrated Parts

Manual Text Revision History

<u>Rev</u>	<u>Dash</u>	<u>Date</u>	<u>Change Description</u>
-1		11/20/02	Engineering Release
A	01	12/09/02	Manufacturing Release

This document contains advanced service information that is company confidential. Due to the proprietary nature of the contents the document should be viewed by GE personnel only. This manual may not be reproduced, in whole or in part, without the written permission of GE OEC Medical Systems, Inc.

OEC is a registered trademark of GE OEC Medical Systems, Inc. Other product and company names mentioned herein are the property of their respective owners.

The contents of this document are accurate at the time of publication. However, changes in design and additional features can, at any time, be incorporated in the hardware and software and may not be reflected in this version of the document. Contact GE OEC Technical Support for clarification, if discrepancies arise.

© 2002  
GE OEC Medical Systems, Inc.  
384 Wright Brothers Drive  
Salt Lake City, Utah 84116  
U.S.A.  
(801) 328-9300

**WARNING**

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRICAL SHOCK, MECHANICAL OR OTHER HAZARDS.

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